VERMEIREN

LiNX 200 controller

INSTRUCTION MANUAL
MODE D'EMPLOI
GEBRUIKSAANWIJZING
GEBRAUCHSANWEISUNG
ISTRUZIONI PER L'USO
MANUAL DE INSTRUCCIONES
INSTRUKCJA OBSŁUGI

EN Contents

1. Use 2

1.1. Turn wheelchair On/Off 3

1.2. Sleep mode 3

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OPERATOR CONTROL LiNX 200



1. Use

The operator control built into your electric wheelchair enables you to control all the driving, steering, braking and operating processes and to control other adjusting motors (lifting column, seat and back adjuster). The wheelchair's electrical unit and electronics are constantly being monitored internally. Any fault in the electronics is indicated on the status indicator (1) (see the chapter on troubleshooting).


Vermeiren is responsible for the changes in the software. For changes in the software contact Vermeiren.

⚠ WARNING: Risk of burns – Be aware that the surface of the operator control can potentially get hot when it is exposed to strong sunlight for long periods.

- 1 = "ON / OFF" button, Status indicator
- 2 = Battery charge indicator
- 3 = Speed switch: left-slower; right-faster
- 4 = "Horn" button
- 5 = Joystick
- 6 = Drive function
- 7 = Adjusting function selector
- 8 = Drive/actuator status
- 9 = Hazard lights and left indicator
- 10 = Driving lights and right indicator
- 11 = Software installation
- 12 = Charger connector

Without lightning

With lightning

Press the "ON/OFF" button (1). The charge indicator (2), which also indicates the batteries' charge level, lights up for a moment. If all the lights come on, the batteries are sufficiently charged. If the lights do not all come on, you must adjust your driving activities to suit the low capacity of the batteries. If the battery charger display lights do not come on, check the plug connections of the electronic system.

Now select the driving program you want by turning the speed switch (3), left turn to go slower and right turn to go faster. Select the drive function with the drive function button (6).

To set the wheelchair in motion, push the joystick (5), forwards slowly. The speed will increase when you push the joystick further forwards. If you want to turn left or right, simply push the joystick in the desired direction. To reverse, pull the joystick backwards from the central position.

⚠ Note that, when pressing the "ON/OFF" button, the joystick stays in the central position for at least five seconds. It is programmed this way for technical and safety reasons, to prevent the joystick being pushed at the same time as the wheelchair is being switched on. If both processes are activated at the same time, the steering unit locks itself and can only become operational again after being switched off completely.

You can sound an acoustic warning signal with the "Horn" (4) button. To activate the adjust functions, press button (7) upwards/downwards until the corresponding function is shown on the drive/actuator status (8).

Press and hold the "Left indicator" (9) or "Right indicator" (10) buttons to signal a change in your driving direction. The driving lights are switched on or off by short pressing the "Lights" (10) button.


The hazard warning lights, used in critical situations, can be switched on or off by short pressing the “Hazard warning lights” **(9)** button.

To brake while driving, simply push the joystick to its central position. Move the joystick slowly for gradual braking. Simply release the lever for a quick stop; the wheelchair will then stop as quickly as possible.

Turn wheelchair On/Off

Turn the power on:


Press the On/Off button. The On/Off button will light green and all LED's of the battery indicator turn on and display the current battery status.

If there is a fault with the system when powering up, the status indicator will indicate the fault with a series of red flashes. If the fault is one that prevents the system from driving, then the battery gauge will flash continuously.

If the wheelchair is turned on while the joystick is out of neutral, the wheelchair will not drive. Release the joystick back to neutral in 5 seconds and the fault will disappear. If the joystick is not released back to neutral (centre) within 5 seconds, the fault becomes a Module fault (the status indicator is flashing red), and the wheelchair will not drive. To clear this fault, loosen the joystick to the neutral position and turn the operator control off and then back on.

Turn the power off:

Press the On/Off button. The system will power down and the status indicator will switch off.

Emergency stop:

If the user needs to stop the wheelchair quickly, or stop a seating motion quickly, the On/Off button can be pressed to perform an EMERGENCY STOP. The wheelchair will come to a halt quickly.

Sleep mode


The operator control will go to sleep after a period of inactivity of a quarter.

When the system sleeps, it is partially turned off to reduce energy consumption.

To reactivate the operator control:

Press any button on the operator control or move the joystick. The system turns on.

Using the battery charge indicator

The battery charge indicator (located above the on/off button) is used to indicate that your wheelchair power is on and provides an estimate of the remaining battery capacity.

Any green LED's indicate that your batteries are well charged.


If only amber and red LED's are light, the batteries are moderately charged. Recharge your batteries before you make a long trip.

If only red LED's are light, the batteries are running out of charge. Recharge as soon as possible.

The following table indicates what the battery charge indicator will display.

Display	Description	Meaning	Remarks
	All LED's dark.	System off.	
	All LED's are lit and the green LED's are flashing.	System on. High voltage warning: battery voltage level has risen above the high voltage warning set-point.	Warning: battery too full. Slow down if you are driving down a slope.
	All LED's are lit.	System on.	LED's show the charge level.: Full battery (between 81 and 100% of full capacity)
	1 Green led is lit.	System on.	Battery level is decreasing.
	Only red and amber LED's lit.	Low battery level.	Battery half full. Start the return journey.
	Red and amber LED's lit.	Very low battery level.	Battery low (between 0 and 20% of full capacity). Recharge soon.
	Only Red LED lit.	Battery almost empty.	Recharge now.
	1 Red Led is flashing.	Battery empty. Low voltage warning: battery voltage level decreased below its low voltage warning set-point.	Recharge immediately.
	First red LED will flash on the battery indicator. The status indicator will display a flash code 2 or 7. The horn will sound once every 10 seconds.	Battery voltage decreases below the battery cut-off voltage.	Battery is fully discharged, and it is good practice to stop discharge at this point. The display and wheelchair shut down automatically.

Battery charging


Plug the battery charger into the remote's XLR socket.

The battery charger indicator will indicate the system is being charged by cycling between a left to right chase sequence, and then displaying the approximate battery charge state at the end of the chase sequence.

Driving is prevented (inhibited) while the system is being charged.


The LiNX 200 system does not have to be powered up when charging the battery, however, if it is not powered up, then the battery gauge will not display the charging state / chase sequence.

The battery charger connection is below the petrol pump symbol.

Adjusting the driving speed


The speed switch allows the user to limit the maximum speed of the wheelchair (that is the speed when the joystick is fully deflected) their preference and environment.

The speed switch offers 10 discrete steps between the lowest speed (switch to the left) and the highest speed (switch to the right).


As a visual reminder, a speed symbol is positioned above the speed switch to indicate the low and high positions.

Using the horn

Press the horn button (located above the joystick). The horn shall sound as long as the button is pressed.

The status indicator

The status indicator is integrated in the ON/OFF button. When the LiNX 200 system is not powered up, the status indicator is not lit.

When the LiNX 200 system is powered up, and there are no faults with the system, the status indicator will be lit green.

If, when powered up, there is a fault with the system, then the status indicator will flash red. The number of flashes will indicate the type of error. For flash codes refer to paragraph 2: Troubleshooting.


Using the Joystick

Moving the joystick will drive the electric wheelchair in that direction. The amount of joystick movement will determine the speed that the powerchair will move in that direction.

Make sure that the joystick is in the neutral (central) position when pressing the on/off button, otherwise the electronics will be locked. To take this block away let go the joystick until he goes back in neutral position.

Drive inhibit indication


Drive inhibit mode is indicated by the battery charger indicator with a right-to-left chase sequence.

The chase sequence starts with the green LED on the right-hand side, and by one-by-one, each LED will switch on and then off. When the sequence completes at the left-most red LED, it begins again at the righthand side.

The chase sequence continues until the error condition has been cleared.

Changing drive function

The wheelchair's drive function can be selected with the drive function select button, which is located below the power button on the left-hand side of the remote module.


Press the top of the button to select the next drive function. Press the bottom of the button to select the previous drive function.

The drive function selected is indicated on the drive function indicator, which is located to the right-hand side of the drive select button.

The indicator has three LED's:

1. when drive function 1 is selected, the bottom indicator will be lit;
2. when drive function 2 is selected, the bottom and middle indicators will be lit;
3. when drive function 3 is selected, all indicators will be lit;

Changing adjusting function

The wheelchair's adjusting functions can be selected with the adjusting function select button, which is the rocker button located on the right-hand side of the operator control.

Press the top of the button to select the next adjusting function. Press the bottom of the button to select the previous adjusting function. Note that only adjusting functions that have been programmed will be available for selection.


The selected adjusting function is displayed on the drive/actuator status indicator as shown below. The display shown only the available/programmed functions in your wheelchair.

Adjust function	Display	Adjust function	Display
Seat inclination		Backrest inclination	
Lift		Leg rest left	
Leg rest right		Leg rests simultaneous	
Unspecified		None	


Drive and adjusting function indications

Drive slow-down


To notify the user that a drive slow-down is active (only for the lift function), the drive wheel LED and the lift function LED's are slowly pulsed on and off. The LED's pulse for the duration of driving or adjusting function demand.

Lights (if buttons available)

There are two lighting control buttons, which are situated below the horn button. The left-hand button controls the left indicator and the hazard lights. The right-hand button controls the right indicator and the driving lights. The operation of these is described below.


Hazard lights

To switch on the hazard lights, press and quickly release the left-hand lighting button once. LED's behind the left-hand button and right-hand button will flash on and off for the duration the hazard lights are operating.

To switch off the hazard lights, press the left-hand or right-hand lighting button.

Indicator lights


To switch on the left indicator, press and hold the left-hand lighting button once. The LED behind the left-hand button will flash on and off for the duration the left indicator lights are operating.

To switch on the right indicator, press and hold the right-hand lighting button once.

The LED below the right-hand button will flash on and off for the duration the right indicator lights are operating.

To switch off either indicator, press the left-hand or right-hand lighting button.

Driving lights


To switch on the driving lights, press and quickly release the right-hand lighting button once. The LED behind the right-hand button will remain on for the duration the driving lights are operating.

To switch off the position lights, press the right-hand or left-hand lighting button.


Using the lock

The lock function is used, primarily, to restrict who can use the system, but also can help prevent unintentional use of the controls for when the system is not required for any length of time. When a system is locked (see below), the system is powered down, and the user controls are not responsive. If the power button is pressed when the system is locked, the locked status is displayed to the user by the battery gauge.

To unlock the system, an unlock sequence must be performed (see below) by the user within a specific time frame. If the sequence is not performed correctly, within the time frame, the system will remain locked and the system will power down again.

To **lock** the system, press and hold the power button for 4 seconds. When entering the locked state, the battery gauge will indicate the transition by flashing LED's 1, 3, and 5 (far left, middle, and far right) 3 times.

To **unlock** the system, press the power button once, and then, press the horn button twice — the horn button must be pressed twice within 10 seconds of pressing the power button.

If the user implements the unlock sequence incorrectly, or the power button is pressed again before the unlock sequence is complete, the system will return to the locked state.


During an unlock attempt, the battery gauge will indicate the system is in a locked state by flashing LED's 1, 3, and 5 (far left, middle, and far right) until either the system is powered off, unlocked, or the Sequence Timeout ± 15 s is reached and the controller switched off automatically.

2. Troubleshooting

The following diagnoses refer to all the electronics. The various modules for steering, power have their own light signals indicating the status of the module concerned (steady light = OK; flashing light = defective). The flashing code is read from the status indicator on the operator control.

If, when powered up, there is an error with the system, then the status indicator will flash red. The number of flashes will indicate the type of error. These is described in the table below.

This fault analysis aids problem analysis and problem solving by indicating the possible faults. If faults arise which cannot be identified and remedied with the aid of this analysis, kindly contact the manufacturer. We wish to point out that in any case of improper handling or if the delivered configuration has been changed, we cannot guarantee the accuracy of the fault indications.

Number	Problem/Fault	Checks
1	Remote / joystick error	<ol style="list-style-type: none"> 1. Check that the joystick is placed in the neutral position when switch on the wheelchair. 2. Check that the joystick is not damaged and function very well.
2	Network or configuration error	<ol style="list-style-type: none"> 1. Check all connection points, cables. 2. Verify software. 3. If you cannot found any fault contact your specialist dealer.
3	Left motor (or its connection) is defective.	<ol style="list-style-type: none"> 1. Check whether the plugs of <u>both</u> motors have been inserted properly. 2. Check the plug contacts of <u>both</u> motors for corrosion or damage. 3. Check <u>both</u> motors. Unplug the motors and measure the plug connection with an Ohm meter. If you obtain readings of more than 1 Ohm or less than 100 milliOhms, the motor is defective. 4. Check the resistance of the motor. Measure every contact using an Ohm meter. <p>⚠ A Fault is also indicated if the motors are disconnected. Connect again and restart the system.</p>
4	Right motor (or its connection) is defective.	As described above.
5	Left parking brake (or its connection) is defective or neutral active.	<ol style="list-style-type: none"> 1. Check whether the motors' plugs have been inserted properly. 2. Check the plugs for corrosion or damage. 3. Check the parking brakes. Measure the resistance of the connections by means of an Ohm meter. If the resistance value is above 100 Ohms or beneath 20 Ohms, the parking brake is probably defective.
6	Right parking brake (or its connection) is defective.	As described above.
7	Module error (other than Remote)	One of the wheelchair parts can be damaged except the operator control itself.

3. Electromagnetic compatibility

The controllers have been tested on typical, representative vehicles to confirm compliance with the following appropriate EMC standards: EN 12184: 2014, ISO7176-21: 2009.

4. Safety instructions

- Do not install, maintain, or operate this equipment before you have read and understood all the instructions and all the manuals for this product and all the other products that you use or install together with this product. Follow the instructions of the manuals. If you do not follow all instructions, injury or damage can be the result.

- Do not try to open or disassemble any case - there are no user-serviceable parts inside.
- The operator has the responsibility to keep the wheelchair in a good safe operating condition. To protect all the components (for example the cables) from damage, the operator must fasten them in optimum positions.
- Do not touch the connector pins. If you touch the pins, they can become dirty or they can be damaged by electrostatic discharge.
- Immediately turn the controller off and consult your service agent if the vehicle:
 - o is damaged;
 - o does not behave the same every time;
 - o does not respond normally, the way you expect it to;
 - o becomes hotter than normal;
 - o smokes;
 - o does not change its speed when you adjust the speed dial;
 - o displays a fault on its fault indicator and the controller does not perform normally.
- Turn the controller off:
 - o when you do not use it;
 - o before you get in or get out of the wheelchair;
 - o if your vehicle driven by itself or against your will. When you turn the controller off the vehicle will halt.
- In the case of an emergency while the vehicle is driving, or a adjusting motion is operating, press the On/Off button to perform an emergency stop and turn the controller off.
- If there is a risk of collision with a person or object in close proximity, use the Joystick and/or speed dial to reduce the speed of the wheelchair.
- Do not drive the vehicle if the controller indicates that the battery is low, since the wheelchair may stop operating and the user may become stranded. If the battery becomes completely empty, the vehicle will stop suddenly and the battery may be damaged.
- Make sure that the battery charger that is used with the vehicle has a drive inhibit function that is correctly connected for use with the controller. This drive inhibit feature is intended to prevent the chair from driving while the battery is being charged. If you are not sure, ask your dealer.
- Do not use battery chargers with a higher current rate that is specified in the manual of your wheelchair.
- If operator of the wheelchair is left with limited or no mobility for any reason (for example, because the wheelchair loses electric power or breaks down), it is important that they can seek assistance from wherever they may be.
- Go downhill slowly. Decrease the speed of the wheelchair when going downhill.
- The controller can cause the vehicle to come to a sudden stop. If this can be dangerous to the operator, the installer must install a seat belt, and the operator must wear this seat belt.
- Operation of a vehicle on steep slopes can be dangerous. Before you drive up or down a slope, make sure that the slope does not exceed the maximum safe slope of your wheelchair.
- Do not use the park brake release on a slope.
- Make sure that the controller does not become colder or hotter than the minimum and maximum temperatures: -25°C to +50°C.
- Most electronic equipment is influenced by Radio Frequency Interference (RFI). Be careful when portable communications equipment is used in the area around such equipment.
- Performance adjustments must only be made by healthcare professionals, or by persons who completely understand the programming parameters, the adjustment process, the configuration of the vehicle, and the capabilities of the driver. Wrong settings can make the vehicle uncontrollable or unstable. An uncontrollable or unstable vehicle can cause an unsafe situation such as a crash, with the risk of serious injury to the driver or bystanders, or damage to the vehicle or surrounding property.
- Performance adjustments must only be made indoors or outdoors in dry conditions.
- If the wheelchair has not been fitted with a hand guard for the remote module, the user should be made aware that their hand will not be protected from crushing, when, for example, maneuvering towards a table.
- Users should be aware that the surface of the controller can potentially get hot when exposed to strong sunlight for long periods.
- The XLR connector on the controller is to be used exclusively for the intended purpose of charging the battery and/or programming the controller. Warranty will be voided if any unauthorized device is connected to this port.

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Vermeiren Suisse S.A.